

Reference level, stability and sample material

Test	Reference level ¹ (adults)		Stability and sample material							Comments
	Conventional	System International	Serum	fresh capillary and fresh venous blood	Heparin Blood	EDTA Blood	Heparin Plasma	EDTA Plasma ⁴		
Bilirubin	< 1.2 mg/dL	< 21 µmol/L	2 h	Immediately	2 h	2 h	2 h	2 h	EDTA and Heparin blood to be kept in the dark. The Reflotron Bilirubin test can not be used on newborn because of special rheological and physiological characteristics of the sample material	
Cholesterol	< 200 mg/dL	< 5.2 mmol/L	4 h ² 12 h ³	Immediately	4 h ² 8 h ³	4 h ² 8 h ³	4 h ² 12 h ³	4 h ² 12 h ³	Shake the sample before performing tests. Do not freeze specimens	
Creatinine	< 1.1 mg/dL < 0.9 mg/dL	< 97 µmol/L < 80 µmol/L	24 h	Immediately	6 h	6 h	24 h	24 h	Urine as sample please see packaging insert	
Glucose	60 - 109 mg/dL	3.3 - 6.05 mmol/L	2 h	Immediately	10 min	10 min	2 h	2 h	Separate serum from the cellular components immediately after coagulation - but no later than 1/2 hour after collection of the blood sample. Perform the glucose determination within 2 hours	
Hemoglobin	♂ < 14 - 17.5 g/dL ♀ < 12.3 - 15.3 g/dL	< 8.7 - 10.9 mmol/L < 7.6 - 9.5 mmol/L	×	Immediately	24 h	24 h	×	×	Heparin and EDTA blood after storage give a good shake	
HDL HDL Cholesterol	< 40 mg/dL or < 1.04 mmol/L low HDL Cholesterol (high risk for CHD) ≥ 60 mg/dL or ≥ 1.56 mmol/L high HDL Cholesterol (high risk for CHD) CHD = Coronary Heart Disease		×	×	×	×	×	after 24 h +4%	Only fresh EDTA plasma	
K+ (Potassium)	Serum Plasma	3.6 - 5.0 mmol/L 3.5 - 4.6 mmol/L	>7 d	×	×	×	>7 d	×	Separate serum from the clot within one hour; if heparinized plasma is used, the anticoagulated blood should be centrifuged and separated from the cells immediately after collection. In closed containers potassium is stable in serum or plasma for at least one week at +4°C to +25°C	
TG Triglycerides	≤ 200 mg/dL	≤ 2.30 mmol/L	8 h ² 24 h ³	Immediately	8 h	8 h	8 h ² 24 h ³	8 h ² 24 h ³	Kept in closed container. Do not freeze specimens	
UA Uric Acid	♂ 3.4 - 7.0 mg/dL ♀ 2.4 - 5.7 mg/dL	200 - 420 µmol/L 140 - 340 µmol/L	3 d ² 7 d ³	Immediately	8 h	×	3 d ² 7 d ³	×	Kept in closed container	
Urea	< 65Y < 50 mg/dL > 65Y < 71 mg/dL	< 8.3 mmol/L < 11.9 mmol/L	7 d	Immediately	8 h	8 h	7 d	7 d	Kept in closed container	
ALP. Alkaline Phosphates	♂ 106 - 265 U/L ♀ < 50Y normal weight 83 - 223 U/L ♀ > 50Y over weight 91 - 258 U/L	1.77 - 4.42 µkat/L 1.39 - 3.72 µkat/L 1.52 - 4.30 µkat/L	3 d 7 d ³	Immediately	8 h	×	1 d	×	Kept in closed container	
AMYL Amylase	< 100 U/L	< 1.65 µkat/L	7 d ² 1 m ³	Immediately	8 h	×	7 d ² 1 m ³	×	Kept in closed container	
P-AM Pancreatic Amylase	< 53 U/L Random Urine 325 U/L or 5.42 µkat/L	< 0.89 µkat/L	7 d	Immediately	8 h	×	7 d	×	Kept in closed container. For urine samples use reference level: Random urine 2d ² ≥ 10d ³	
CK	♂ 24 - 195 U/L ♀ 24 - 170 U/L	0.4 - 3.25 µkat/L 0.4 - 2.83 µkat/L	24 h ² 7 d ³	Immediately	8 h	×	24 h ² 7 d ³	×	Kept in closed container	
GGT	♂ 10 - 71 U/L ♀ 6 - 42 U/L	0.17 - 1.19 µkat/L 0.10 - 0.70 µkat/L	7 d	Immediately	8 h	8 h	7 d	7 d	Kept in closed container	
GOT (AST)	♂ < 40 U/L ♀ < 33 U/L	< 0.67 µkat/L < 0.55 µkat/L	4 d ² 7 d ³	Immediately	1 h	×	4 d ² 7 d ³	×	Kept in closed container	
GPT (ALT)	♂ < 41 U/L ♀ < 32 U/L	< 0.68 µkat/L < 0.53 µkat/L		Immediately	1 h	×	3 d ² 7 d ³	×	Kept in closed container	

¹ testing Temperature 37°C

² kept in closed container at +20 - 25°C

³ kept in closed container at +4 - 8°C

⁴ plasma separation to be done immediately



Life needs answers